

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330021

Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

Adult (M) cotton, velcro, grey: 330011

Adult (M) cotton, velcro, red : 330013

Adult (M) navy, with straps: 330112

Adult (M) red, with straps: 330113

Adult (L) cotton, velcro, navy: 330002

Child/small adult (S), cotton, velcro, grey: 330021

Features

CE Marking	NONE	Medical device	No
Bar code	3700446006822	HS Code	9018901000
Warranty (months)	6	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Packed unit / Height (cm)	3
Packed unit / Length (cm)	15	Packed unit / Weight	0.138

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 30/04/2026 00:44:08 – For information purposes only, can be subject to change