

VAQUEZ-LAUBRY CLASSIC VELCRO CUFF WITH BLADDER/DARK BLUE/LARGE ADULT



Ref. 330112

Description

Made of polycotton with inner bladder and dual tubing (10 cm & 60 cm) - also suited for Vaquez-Laubry® with full-grain leather manometer holder, and Velcro fastening.

Adult (M) cotton, velcro, grey: 330011

Adult (M) cotton, velcro, red : 330013

Adult (M) navy, with straps: 330112

Adult (M) red, with straps: 330113

Adult (L) cotton, velcro, navy: 330002

Child/small adult (S), cotton, velcro, grey: 330021

Features

| | | | |
|-------------------|---------------|----------------|------------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3700446031886 | HS Code | 9018901000 |
| Warranty (months) | 6 | Country | France |

DIMENSIONS AND WEIGHT

| | | | |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.225 |

SPECIFIC STORAGE CONDITIONS

| | | | |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |

OTHER FEATURES

| | |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

📅 30/04/2026 01:16:51 – For information purposes only, can be subject to change