

PEDIATRIC CUFF 18-26 CM FOR 4VS MONITOR

Ref. 551621



Description

Pediatric cuff for PNI module of SPENGLER 4VS vital signs monitor.
Dimensions (18-26cm)

30/04/2026 00:47:26 – For information purposes only, can be subject to change

PEDIATRIC CUFF 18-26 CM FOR 4VS MONITOR



Ref. 551621

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446033057	HS Code	9018901000
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.049
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/04/2026 00:47:26 – For information purposes only, can be subject to change