

CATARACT SURGERY SET

Ref. IBOI403



30/04/2026 00:25:51 – For information purposes only, can be subject to change

CATARACT SURGERY SET

Ref. IBOI403



Description

This box contains:

- Stainless steel box, 20 x 10 x 5 cm -IS
- Bowman blepharostat, 7 cm
- Iridectomy scissors, curved, 10 cm -IS, Holtex
- *Vannas corneal scissors, 10 cm, straight, foam-coated
- Castroviejo scissors, corneal, curved, soft, 10 cm
- Graefe cataract knife 12 cm/20 mm
- Graefe iris hook, fine
- Micro Halstead forceps, 10 cm, straight, A/G, Holtex
- *Barraquer fixation forceps, 1 x 2 oblique claws
- Elschmig Fixing Forceps, A/G, oblique
- *Bishop-Hartmann Dissection Forceps, 8 cm, straight, S/G
- *Bishop-Hartmann Dissection Forceps, 8 cm, straight, A/G
- Graefe Dissection Forceps, fine, iris, 10 cm, curved, S/G
- *Wecker Hemostatic Forceps for Tensioning Sutures, 5 cm, straight
- *Castroviejo Needle Holder, foam, 14 cm, straight
- Castroviejo Iris Spatula
- De Wecker Iris Spatula
- Silicone Mat, 20 x 30 cm

Features

CE Marking	NONE	Medical device	No
Bar code	3661040082660	HS Code	90189084
Warranty (months)	24	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Packed unit / Height (cm)	5
Packed unit / Length (cm)	21	Packed unit / Weight	1
Bare Unit / Weight (kg)	1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

30/04/2026 00:25:51 – For information purposes only, can be subject to change