

# LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

Ref. [MEEMER03\\_B](#)

**EMERGENCY'S™**



GARANTIE  
2 ANS  
2 YEARS  
GUARANTEE



30/04/2026 00:16:13 – For information purposes only, can be subject to change

## LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

**EMERGENCY'S™**

Ref. [MEEMER03\\_B](#)

**Intended to transport the most fragile + relief materials**

### Description

- Four compartments for interior supports and high-capacity compartment.
- Transparent pouch in the inner flap, designed to see and organize the material easily, thanks to the rubber-elastic that guarantees a correct fixation- Five outer pockets, two side pockets, one on the back and one on the top and a front pocket that can be divided internally thanks to rubber pockets and bands
- Base of the waterproof bag
- High resistance zips
- Reflective strips- Reinforced handles and shoulder straps
- External card holder
- Dimensions: 55 x 31.5 x 33 cm, capacity 38 L
- Internal dimensions: 40 x 20 x 24 cm
- Material 600D Polyester
- Weight 1.76 kg
- Sold empty
  
- Available in : Red (MEEMER03\_R), Blue (MEEMER03\_B) and Orange (MEEMER03\_O)

 30/04/2026 00:16:13 – For information purposes only, can be subject to change

# LARGE CAPACITY EMERGENCY BAG (38L / 55X31.5X24CM) - BLUE

**EMERGENCY'S™**

Ref. [MEEMER03\\_B](#)

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564542588	HS Code	420299000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	28	Bare Unit / Width (cm)	31.5
Packed unit / Height (cm)	20	Bare Unit / Height (cm)	33
Packed unit / Length (cm)	50	Bare Unit / Length (cm)	55
Packed unit / Weight	1.8	Bare Unit / Weight (kg)	1.76
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Blue	Packed unit / Quantity	1

30/04/2026 00:16:13 – For information purposes only, can be subject to change