

REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011G



CE



📅 02/05/2026 11:35:00 – For information purposes only, can be subject to change

REUSABLE NIBP FOR VITALIMATE MINI

Ref. MONI011G

Vitalimate MINI reusable NIBP cuffs

Description

The Vitalimate MINI NIBP cuffs are comfortable and come in all sizes for all types of patients.

- Reusable adult NIBP cuff (27-35cm) : ref MONI011
- Reusable large adult NIBP cuff (34-43cm) : ref MONI011A
- Small reusable adult NIBP cuff (20.5-28 cm): ref MONI011G
- Reusable pediatric NIBP cuff (14-21.5 cm): ref MONI011C
- Reusable neonatal NIBP cuff (10-15 cm): ref MONI011D
- Reusable adult NIBP thigh cuff (42-54 cm): ref MONI011B

Disposable neonatal cuffs are also available for infants.

Features

CE Marking	CE	Medical device	Class I
Bar code	6948309002575	HS Code	9018193090
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085
Bare Unit / Weight (kg)	0.08		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	10	Maximum humidity	93

OTHER FEATURES

📅 02/05/2026 11:35:00 – For information purposes only, can be subject to change