

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC

ELITE BAGS



04/05/2026 23:42:07 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC



Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner

Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening

- Material: Polyester gray bitone
- Weight: 0.12 kg.

- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03_N), blue jeans (MEDIAB03_B) or pink (MEDIAB03_RC)

- Sold empty

 04/05/2026 23:42:07 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC



Features

CE Marking	NONE	Medical device	No
Bar code	8436564542465	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	17
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.012

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	30	Maximum humidity	85

OTHER FEATURES

Colors	Pink	Packed unit / Quantity	1
--------	------	------------------------	---

📅 04/05/2026 23:42:07 – For information purposes only, can be subject to change