

Ref. 332012



CE



04/05/2026 23:17:35 – For information purposes only, can be subject to change

Ref. 332012


Description

Cotton cuff with highly resistant rubber bladder.
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
 Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028480	HS Code	9018901000
Warranty (months)	6	Country	Indonesia

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.21
Bare Unit / Weight (kg)	0.21		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--

 04/05/2026 23:17:35 – For information purposes only, can be subject to change