



LED FIBER OPTIC OTOSCOPE - SMARTLED 5500

Ref. 570517B



04/05/2026 23:16:24 – For information purposes only, can be subject to change

LED FIBER OPTIC OTOSCOPE - SMARTLED 5500

Ref. 570517B



The SMARTLED 5500 otoscope has been developed with a lens to allow exceptional diagnosis.

Description

LED lighting with Optical fibreThe Spengler otoscope SMARTLED 5500® is using LED lighting with optical fibre allowing a direct, homogeneous illumination and without reflections of the auricular canal and the eardrum. Luminance 85K LUX Bimaterial head in ABS and Stainless steel, screwed head for a safe fixation.Ergonomic ABS grip. End clamp with built-in switch button.Optical fibre for optimal transmission of the light beam. x 3 magnification.Supplied in a case with two LR6 AA batteries, 8 speculums (two of 2.5 mm and four of 4 mm).

Available in 3 colors: black (570517), coral (570517C), blueberry (570517B) and also in a paediatric version (light green 570517P)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446041717	HS Code	90189042
Warranty (months)	36	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	4
Packed unit / Length (cm)	20	Packed unit / Weight	0.385
Bare Unit / Weight (kg)	0.075		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	Blueberry	Packed unit / Quantity	1

📅 04/05/2026 23:16:24 – For information purposes only, can be subject to change