

IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_RC




CE



05/05/2026 00:00:22 – For information purposes only, can be subject to change

IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_RC



Description

- Stethoscope with double metal horn adult (Ø: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073_B), Black (DH7073_N), Red (DH7073_R) or Pink (DH7073_RC) for the paediatric version (Ø: 35/30 mm)

Features

| | | | |
|------------------------------------|---------------|---------------------------|------------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661041000106 | HS Code | 9018908400 |
| Warranty (months) | 12 | Country | China |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 12 | Packed unit / Height (cm) | 4 |
| Packed unit / Length (cm) | 25 | Packed unit / Weight | 0.15 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 70 |
| Minimum humidity | 20 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Colors | Pink | Packed unit / Quantity | 1 |

📅 05/05/2026 00:00:22 – For information purposes only, can be subject to change