

IDEAL+ DOUBLE HEAD ADULT, BLUE


Ref. DH7073_R



05/05/2026 00:06:23 – For information purposes only, can be subject to change

IDEAL+ DOUBLE HEAD ADULT, BLUE

Ref. DH7073_R



Description

- Stethoscope with double metal horn adult (Ø: 45/30 mm)•
- Anti-cold ring and horn of the color of the tube•
- Aluminum lyre incorporated into the tubing
- Also available in single pavilion for adults

Available in Blue (DH7073_B), Black (DH7073_N), Red (DH7073_R) or Pink (DH7073_RC) for the paediatric version (Ø: 35/30 mm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3661041000090	HS Code	9018908400
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.35
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Colors	Red	Packed unit / Quantity	1

📅 05/05/2026 00:06:23 – For information purposes only, can be subject to change