

# EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Ref. [MEEMER02\\_N](#)



GARANTIE  
2 ANS  
2  
2 YEARS  
GUARANTEE



📅 05/05/2026 00:12:50 – For information purposes only, can be subject to change

## **EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)**

**EMERGENCY'S™**

Ref. MEEMER02\_N

**Practical and lightweight first aid bag with two adjustable internal dividers and a zippered pocket in the flap to facilitate the organization of sanitary equipment**

### **Description**

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High resistance zip closures

-

- Reflective strips- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material 600D Polyester
- Weight 1.29 kg

-

Available in waterproof version

- Sold empty

- Available in blue (MEEMER02\_B), yellow (MEEMER02\_J), black (MEEMER02\_N), orange (MEEMER02\_O), pink (MEEMER02\_RC) and red (MEEMER02\_R).

# EMERGENCY EMERGENCY BAG LIGHT (22 L / 44X25X27CM) - BLUE (EM13.014)

**EMERGENCY'S™**

Ref. [MEEMER02\\_N](#)

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564542977	HS Code	420299000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	24	Bare Unit / Width (cm)	26
Packed unit / Height (cm)	22	Bare Unit / Height (cm)	19
Packed unit / Length (cm)	30	Bare Unit / Length (cm)	31.6
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.29
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

📅 05/05/2026 00:12:50 – For information purposes only, can be subject to change