

**MEDIPLUGS OPERATING ROOM SHOES,
NON-PERFORATED, 35, NAVY BLUE**

Ref. WSAM_39V

medi **Plogs**



📅 05/05/2026 00:40:56 – For information purposes only, can be subject to change

MEDI PLOGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM_39V



Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

 05/05/2026 00:40:56 – For information purposes only, can be subject to change

**MEDIPLUGS OPERATING ROOM SHOES,
NON-PERFORATED, 35, NAVY BLUE**



Ref. WSAM_39V

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | NONE | Medical device | No |
| Bar code | 3661040099972 | HS Code | 64029993 |
| Warranty (months) | 12 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 1.16 |
| Bare Unit / Weight (kg) | 1 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Green | Shoe Size | 39 |

📅 05/05/2026 00:40:56 – For information purposes only, can be subject to change