

**MEDIPLUGS OPERATING ROOM SHOES,  
NON-PERFORATED, 35, NAVY BLUE**

Ref. WSAM\_43L


medi **Plogs**



📅 04/05/2026 23:37:48 – For information purposes only, can be subject to change

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM\_43L



### Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

 04/05/2026 23:37:48 – For information purposes only, can be subject to change

**MEDIPLUGS OPERATING ROOM SHOES,  
NON-PERFORATED, 35, NAVY BLUE**



Ref. WSAM\_43L

**Features**

CE Marking	NONE	Medical device	No
Bar code	3661040100197	HS Code	64029993
Warranty (months)	12	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1.16
Bare Unit / Weight (kg)	1		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	fuchsia	Shoe Size	43

📅 04/05/2026 23:37:48 – For information purposes only, can be subject to change