

DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER

Ref. L20747



06/05/2026 22:44:15 – For information purposes only, can be subject to change

DISPOSABLE SPECULUM DIA. 2.5 MM, BOX OF 250, BLACK, CHILD, SPENGLER



Ref. L20747

Auricular devices for accurate and safe otoscopy.

Description

Box of 250.
 For children, single use,
 Black,
 Ø 2.5 mm
 Universal, compatible with all brands
 Box with dispensing spout for easy access to speculums

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446031268	HS Code	3926909760
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Packed unit / Height (cm)	11
Packed unit / Length (cm)	11	Packed unit / Weight	0.14
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

06/05/2026 22:44:15 – For information purposes only, can be subject to change