

# SPENGLER X LUXAMED EDITION - SMARTLED OTOSCOPE - BLACK

Ref. A1.416.112



14/05/2026 22:38:28 – For information purposes only, can be subject to change

# SPENGLER X LUXAMED EDITION - SMARTLED OTOSCOPE - BLACK

Ref. A1.416.112



## SMARTLED SPENGLER X LUXAMED EDITION® OTOSCOPE

An innovative design otoscope made in Germany, tailored to the needs of healthcare professionals.

### Description

Key features :

Design and sturdiness:

- Lightweight (100gr), ideal for fatigue-free everyday use.
- Aluminium and fibreglass-reinforced handle
- Iconic Spengler colour

Bright illumination :

- Provides 10,000 Lux of brightness for clear observation.
- Colour temperature: 4,250 K

Technological innovation:

- Ring of 6 concentric LEDs positioned at the front for a clear, shadow-free field of vision
- Extended service life with constant current, approximately 100,000 hours.
- Adjustable brightness with a long press on the ON/OFF button
- Automatic switch-off after 3 minutes

Quality and innovation MADE IN GERMANY

Also available in Blueberry, Bronze and paediatric version (green)

Includes: 2 AAA/LR03 batteries + 20 disposable specula (10 in 2.5 mm and 10 in 4 mm)

 14/05/2026 22:38:28 – For information purposes only, can be subject to change

# SPENGLER X LUXAMED EDITION - SMARTLED OTOSCOPE - BLACK



Ref. A1.416.112

## Features

CE Marking	CE	Medical device	Class I
Bar code	4260381965305	HS Code	90189084
Warranty (months)	36	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			

📅 14/05/2026 22:38:28 – For information purposes only, can be subject to change