


SPENGLER X LUXAMED EDITION SMARTLED OTOSCOPE - PEDIATRIC



Ref. A1.516.514



CE



14/05/2026 22:30:05 – For information purposes only, can be subject to change

SPENGLER X LUXAMED EDITION SMARTLED OTOSCOPE - PEDIATRIC

Ref. A1.516.514



SMARTLED SPENGLER X LUXAMED EDITION® OTOSCOPE

An innovative design otoscope made in Germany, tailored to the needs of healthcare professionals.

Description

Key features :

Design and sturdiness:

- Lightweight (100gr), ideal for fatigue-free everyday use.
- Aluminium and fibreglass-reinforced handle
- Iconic Spengler colour

Bright illumination :

- Provides 10,000 Lux of brightness for clear observation.
- Colour temperature: 4,250 K

Technological innovation:

- Ring of 6 concentric LEDs positioned at the front for a clear, shadow-free field of vision

 14/05/2026 22:30:05 – For information purposes only, can be subject to change

SPENGLER X LUXAMED EDITION SMARTLED OTOSCOPE - PEDIATRIC



Ref. A1.516.514

- Extended service life with constant current, approximately 100,000 hours.
- Adjustable brightness with a long press on the ON/OFF button
- Automatic switch-off after 3 minutes

Quality and innovation MADE IN GERMANY

Also available in Blueberry, Black and Bronze.

Includes: 2 AAA/LR03 batteries + 20 disposable specula (10 in 2.5 mm and 10 in 4 mm)

SPENGLER X LUXAMED EDITION SMARTLED OTOSCOPE - PEDIATRIC



Ref. A1.516.514

Features

CE Marking	CE	Medical device	Class I
Bar code	4260381965299	HS Code	90189084
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	11	Packed unit / Height (cm)	4
Packed unit / Length (cm)	24	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.12		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			

📅 14/05/2026 22:30:05 – For information purposes only, can be subject to change