


LITESTICK PENLIGHT SPENGLER

Ref. 570622



CE



 18/05/2026 00:24:43 – For information purposes only, can be subject to change

LITESTICK PENLIGHT SPENGLER

Ref. 570622



Pen-shaped medical light, designed for close-up clinical examinations.

Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009878	HS Code	9018908400
Warranty (months)	24		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.015		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
OTHER FEATURES			
Colors	Green	Packed unit / Quantity	1

📅 18/05/2026 00:24:43 – For information purposes only, can be subject to change