

## LITESTICK PENLIGHT SPENGLER

Ref. 570623



CE



 18/05/2026 00:33:58 – For information purposes only, can be subject to change

# LITESTICK PENLIGHT SPENGLER

Ref. 570623



Pen-shaped medical light, designed for close-up clinical examinations.

## Description

This penlight is dedicated to general auscultation

- LED lighting for cold , intense and white light close to natural light
- Robust, made in aluminium
- light and compact
- on/off button
- 2 LR03 AAA batteries included

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009847	HS Code	9018908400
Warranty (months)	24		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	5	Packed unit / Height (cm)	5
Packed unit / Length (cm)	20	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.015		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
<b>OTHER FEATURES</b>			
Colors	Red	Packed unit / Quantity	1

📅 18/05/2026 00:33:58 – For information purposes only, can be subject to change