

COOL BAG FOR DIABETIQUE PM DIA, SILVER

Ref. MEDIA02_A

ELITE BAGS



19/05/2026 00:08:05 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. [MEDIA02_A](#)

Kit designed with the spectrum of patients with diabetes in mind

Description

- Double insulated compartment to keep products cool
- Includes reusable cold gel
- Material: 420D fabric
- Weight: 130 g
- Exterior dimensions: 19 x 10 x 5 cm / 0.95 l
- Max. recommended capacity: 600 g
- Available in blue (ref: MEDIA02_B) or silver (ref: MEDIA02_A)
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	8436564540744	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	20	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	14	Bare Unit / Length (cm)	19
Packed unit / Weight	0.25	Bare Unit / Weight (kg)	0.13

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Argenté	Packed unit / Quantity	1
--------	---------	------------------------	---

📅 19/05/2026 00:08:05 – For information purposes only, can be subject to change