

# COOL BAG FOR DIABETIQUE PM DIA, SILVER

Ref. MEDIA02\_B

**ELITE BAGS**



19/05/2026 00:02:50 – For information purposes only, can be subject to change

# COOL BAG FOR DIABETIQUE PM DIA, SILVER



Ref. [MEDIA02\\_B](#)

**Kit designed with the spectrum of patients with diabetes in mind**

## Description

- Double insulated compartment to keep products cool
- Includes reusable cold gel
- Material: 420D fabric
- Weight: 130 g
- Exterior dimensions: 19 x 10 x 5 cm / 0.95 l
- Max. recommended capacity: 600 g
- Available in blue (ref: MEDIA02\_B) or silver (ref: MEDIA02\_A)
- Sold empty

## Features

CE Marking	NONE	Medical device	No
Bar code	8436564540690	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	12	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	19
Packed unit / Weight	0.214	Bare Unit / Weight (kg)	0.13

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

Colors	Blue	Packed unit / Quantity	1
--------	------	------------------------	---

📅 19/05/2026 00:02:50 – For information purposes only, can be subject to change