

**Ref. 522711**



20/05/2026 01:08:47 – For information purposes only, can be subject to change

**Ref. 522711**

## Description

Cuff available in :

Adult (M), nylon, velcro, black (522711)

Adult (L), nylon, black (522712)

Child/small adult (S), nylon, black (522713)

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028176	HS Code	9018199560
Warranty (months)	60	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.092
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 20/05/2026 01:08:47 – For information purposes only, can be subject to change