

# MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507210



22/05/2026 03:50:57 – For information purposes only, can be subject to change

# MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507210



## Description

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

Available in :

Blue : Ref 507214

Black: Ref 507210

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009120	HS Code	90181990
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.375
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

📅 22/05/2026 03:50:57 – For information purposes only, can be subject to change