


MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507214



22/05/2026 02:36:30 – For information purposes only, can be subject to change

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507214



Description

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

Available in :

Blue : Ref 507214

Black: Ref 507210

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009137	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.375
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Azure	Packed unit / Quantity	1

📅 22/05/2026 02:36:30 – For information purposes only, can be subject to change