

**DIAGNOSTIC PENLIGHT - WITH STANDARD  
BULB - LUXAMED**

Ref. D1.212.112

**LUXAMED®**



📅 22/05/2026 03:31:14 – For information purposes only, can be subject to change

# DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**

Ref. D1.212.112

Penlight for eye and throat diagnostic., long life technology.

## Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

- White: D1.212.912
- Blue: D1.212.212
- Grey: D1.212.312
- Yellow: D1.212.512
- Kids : D1.312.912
- Black: D D1.212.112
- Red : D1.212.412
- Green : D1.212.612

# DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**

Ref. D1.212.112

## Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961741	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.04		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	Black		

📅 22/05/2026 03:31:14 – For information purposes only, can be subject to change