

MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02_B

ELITE BAGS



22/05/2026 02:44:58 – For information purposes only, can be subject to change

MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02_B



Description

Sanitary sports case with practical and comfortable design. Its internal organization totally adaptable thanks to the extractable velcro dividers, allows the rapid location of the equipment in case of medical or sports emergency.

- Quilted handle with velcro closure- Card holder
- Reinforced structure with ribbon
- Isothermal amplestick with 41 elastic strips- Document holder
- Large capacity external pocket with elastic bands and spacious compartments
- Adjustable honeycomb structure interior
- Flat pocket for the document holder
- 2 flat net pockets for instruments
- External pocket for personal oblets
- Measures: 50 x 25 x 31 cm
- Capacity: 33.6L
- Weight: 2.23 kg
- Material: 1200D polyester. Washable
- Recommended max. capacity: 12 kg
- Optional trolley
- Sold empty
- Available in blue (MEMEDI02_B) and black (MEMEDI02_N)

 22/05/2026 02:44:58 – For information purposes only, can be subject to change

MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02_B



Features

CE Marking	NONE	Medical device	No
Bar code	8436564540423	HS Code	4202290000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	31
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	50
Packed unit / Weight	3.55	Bare Unit / Weight (kg)	2.23
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Blue	Packed unit / Quantity	1

22/05/2026 02:44:58 – For information purposes only, can be subject to change