

STETHOSCOPE PULSE II SINGLE

Ref. 507701



23/05/2026 19:08:49 – For information purposes only, can be subject to change

STETHOSCOPE PULSE II SINGLE

Ref. 507701



The essential stethoscope. Minimum weight, maximum performance.

Description

HIGH ACOUSTIC QUALITY

- Single tube with high-performance conduit
- High resolution binaural
- High resolution diaphragm

COMFORT

- Anti-cold ring
- Silicone ear tips
- New headphone design more adapted to the anatomy
- Aluminium earhook to reduce weight

Available in Carbon (507702), Blueberry (507700), Coral (507701) or Powder Pink (507703).

Available with double bell (adult bell and bell).

STETHOSCOPE PULSE II SINGLE

Ref. 507701



Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029876	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.226
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Coral	Packed unit / Quantity	1

📅 23/05/2026 19:08:49 – For information purposes only, can be subject to change