

STETHOSCOPE PULSE II SINGLE

Ref. 507702



23/05/2026 19:31:16 – For information purposes only, can be subject to change

STETHOSCOPE PULSE II SINGLE

Ref. 507702



The essential stethoscope. Minimum weight, maximum performance.

Description

HIGH ACOUSTIC QUALITY

- Single tube with high-performance conduit
- High resolution binaural
- High resolution diaphragm

COMFORT

- Anti-cold ring
- Silicone ear tips
- New headphone design more adapted to the anatomy
- Aluminium earhook to reduce weight

Available in Carbon (507702), Blueberry (507700), Coral (507701) or Powder Pink (507703).

Available with double bell (adult bell and bell).

STETHOSCOPE PULSE II SINGLE

Ref. 507702



Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700446029883 | HS Code | 90181990 |
| Warranty (months) | 60 | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 0.226 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Colors | Carbon | Packed unit / Quantity | 1 |

📅 23/05/2026 19:31:16 – For information purposes only, can be subject to change