

## STETHOSCOPE PULSE II SINGLE

Ref. 507703



23/05/2026 18:02:10 – For information purposes only, can be subject to change

## STETHOSCOPE PULSE II SINGLE

Ref. 507703



The essential stethoscope. Minimum weight, maximum performance.

### Description

#### HIGH ACOUSTIC QUALITY

- Single tube with high-performance conduit
- High resolution binaural
- High resolution diaphragm

#### COMFORT

- Anti-cold ring
- Silicone ear tips
- New headphone design more adapted to the anatomy
- Aluminium earhook to reduce weight

Available in Carbon (507702), Blueberry (507700), Coral (507701) or Powder Pink (507703).

Available with double bell (adult bell and bell).

# STETHOSCOPE PULSE II SINGLE

Ref. 507703



## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446029890	HS Code	90181990
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.226
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Powder pink	Packed unit / Quantity	1

 23/05/2026 18:02:10 – For information purposes only, can be subject to change