

**Ref. 335111**



## Description

Cuffs with dual tube bladder & coupling (60 cm tube with Classic® bulb and valve, and 10 cm with male screw-on metal fitting)

Available in :

Adult (M), nylon, velcro, black (335111)

Adult (L), nylon, black (335101)

25/05/2026 23:47:43 – For information purposes only, can be subject to change

**Ref. 335111**



## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028435	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.245
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

25/05/2026 23:47:43 – For information purposes only, can be subject to change