

Ref. C1.616.914

LUXAMED®



Technologie CCT
2 couleurs de lumière



27/05/2026 01:39:50 – For information purposes only, can be subject to change

=====

Ref. C1.616.914

LUXAMED®

 27/05/2026 01:39:50 – For information purposes only, can be subject to change

Ref. C1.616.914

Innovative dermatoscope in pocket size

Description

Innovative illumination: LED-Ring with 2x5 LEDs in concentric structure for a highly detailed image, homogenous light, switchable with one button press : warm white light or cold white light

- Removable contact plate for easy cleaning
- Maintenance-free, because no lamp replacement necessary
- Stepless zoom area with 10 x magnification for a distortion-free image
- Optimized design of head
- Automatic switchoff after 3 minutes
- Life expectancy of LEDs: approx. 100.000 hours due to ideal operation by constant current drive
- Intensity of illumination: 10.000 Lux
- Colour temperature: ring 1 approx. 3.000 K, ring 2 approx. 6.000 K
- Two-component-handle (aluminium/fibre-glass reinforced),
- Also combinable with LuxaScope Auris CCT LED 2.5 V otoscope-head
- Electronic ON/OFF push-button-switch
- Power supply: 2 x AAA Alkaline batteries, inclusive

Available in :

- Black : C1.616.114
- White : C1.616.914

Ref. C1.616.914

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963837	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	White		

27/05/2026 01:39:50 – For information purposes only, can be subject to change