

-----  
**Ref. E1.423.012**

**LUXAMED®**



📅 27/05/2026 00:36:39 – For information purposes only, can be subject to change

**Ref. E1.423.012**

## Description

LUXAMED® Miller laryngoscope blades with integrated fiber optics are available in several sizes.

	Length (mm)	Width (mm)
Size 00	65	11
Size 0	77	11
Size 1	102	11
Size 2	155	14
Size 3	195	15
Size 4	210	18

- Blade for integrated fiber optics
- Made of high-quality stainless steel for a long service life
- Easy to clean and sterilize as there are no crevices to trap dust
- Maintenance-free
- 5,500 fibers for maximum service life and light intensity
- Optimal illumination with the largest possible cross-section (4 mm diameter)
- Illuminance of approximately 10,500 lux for 3.7V and 2.5V LEDs
- Compatible with all handles compliant with ISO 7376 (eco-friendly standard)
- Autoclavable up to 134°C (approximately 4,000 times)
- Data matrix with all important information for traceability

**Ref. E1.423.012**

**LUXAMED®**

## Features

CE Marking	CE	Medical device	Class I
Bar code	04260381960607	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Size	Size 3		

27/05/2026 00:36:39 – For information purposes only, can be subject to change