
Ref. E1.429.012

LUXAMED®



📅 27/05/2026 00:38:42 – For information purposes only, can be subject to change

Ref. E1.429.012

Description

LUXAMED® Miller laryngoscope blades with integrated fiber optics are available in several sizes.

| | Length (mm) | Width (mm) |
|---------|-------------|------------|
| Size 00 | 65 | 11 |
| Size 0 | 77 | 11 |
| Size 1 | 102 | 11 |
| Size 2 | 155 | 14 |
| Size 3 | 195 | 15 |
| Size 4 | 210 | 18 |

- Blade for integrated fiber optics
- Made of high-quality stainless steel for a long service life
- Easy to clean and sterilize as there are no crevices to trap dust
- Maintenance-free
- 5,500 fibers for maximum service life and light intensity
- Optimal illumination with the largest possible cross-section (4 mm diameter)
- Illuminance of approximately 10,500 lux for 3.7V and 2.5V LEDs
- Compatible with all handles compliant with ISO 7376 (eco-friendly standard)
- Autoclavable up to 134°C (approximately 4,000 times)
- Data matrix with all important information for traceability

Ref. E1.429.012

LUXAMED®

Features

| | | | |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 04260381960560 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Germany |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | 5 | Max temperature | 40 |
| Minimum humidity | 10 | Maximum humidity | 90 |
| OTHER FEATURES | | | |
| Size | Size 00 | | |

 27/05/2026 00:38:42 – For information purposes only, can be subject to change