

# MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)

Ref. 522708



CE  
0123

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
**2**  
2 ANS  
2 YEARS  
GUARANTEE

27/05/2026 23:58:19 – For information purposes only, can be subject to change

## MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)



Ref. 522708

**SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.**

### Description

#### HIGH MEASUREMENT ACCURACY

- Dual hoses for linear pressure build-up

#### COMFORT

- Ambidextrous
- Remote-readable display
- 2.5-meter coiled hose
- Adjustable from 95 cm to 135 cm
- Ergonomic, large-volume bulb for rapid inflation
- Storage basket behind the dial (available as an option)
- Comes with a wheeled stand adjustable from 95 cm to 135 cm

#### CUFF

- Comes with a set of 3 cuffs (S, M, L), nylon

# MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)



Ref. 522708

## Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446031169	HS Code	90189010
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.36
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 27/05/2026 23:58:19 – For information purposes only, can be subject to change