

**POLYESTER AUTOCLAVABLE CUFF
STRAIGHT WITH 1 BLADDER FOR DDM
PNEUMATIC TOURNIQUET**

Ref. GBS101


DDM 



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2
2 ANS
2 YEARS
GUARANTEE



📅 28/05/2026 00:01:13 – For information purposes only, can be subject to change

POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GBS101

Those cuffs are autoclavable!

Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff infant (GBS101)	19 x 4 cm
Cuff child (GBS102)	26.5 x 5 cm
Cuff small adult (GBS101)	35.5 x 6 cm
Cuff adult (GBS103)	53 x 6 cm

**POLYESTER AUTOCLAVABLE CUFF
STRAIGHT WITH 1 BLADDER FOR DDM
PNEUMATIC TOURNIQUET**



Ref. GBS101

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203476	HS Code	90189084
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

28/05/2026 00:01:13 – For information purposes only, can be subject to change