


# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS103



📅 27/05/2026 22:33:51 – For information purposes only, can be subject to change

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS103

Those cuffs are autoclavable!

## Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

Cuff infant (GBS101)	19 x 4 cm
Cuff child (GBS102)	26.5 x 5 cm
Cuff small adult (GBS101)	35.5 x 6 cm
Cuff adult (GBS103)	53 x 6 cm

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS103

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203490	HS Code	90189084
Warranty (months)	24	Country	France

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

### OTHER FEATURES

📅 27/05/2026 22:33:51 – For information purposes only, can be subject to change