

## MULTY'S ELITE BAGS, BLACK

Ref. MEMULT01\_N

**ELITE BAGS**



31/05/2026 18:34:29 – For information purposes only, can be subject to change

## MULTY'S ELITE BAGS, BLACK

Ref. MEMULT01\_N



### Multi-purpose case

## Description

It is designed to offer the professional the necessary support in first aid care in various professions and / or situations, such as: Physiotherapy, First Aid, Sports Activities ... - 1 padded main space, large capacity

with extractable separator- 2

padded pockets of medium capacity

- 1 flat pocket

- 2 net pockets

- 12 medium format elastics

- 8 large format elastics- Handle and shoulder strap- Color: black, blue (yellow shoulder strap)

or gray (pink shoulder strap)

- Weight: 580 gr

- Outdoor measures: 31.5 x 20 x 20 cm / 14 L

- Recommended maximum capacity: 3.5 kg

- Sold empty

- Available in black (MEMULT01\_N) and grey (MEMULT01\_G)

# MULTY'S ELITE BAGS, BLACK

Ref. MEMULT01\_N



## Features

CE Marking	NONE	Medical device	No
Bar code	8436564540409	HS Code	4202290000
Warranty (months)	24	Country	China
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	21	Bare Unit / Width (cm)	20
Packed unit / Height (cm)	18	Bare Unit / Height (cm)	20
Packed unit / Length (cm)	44	Bare Unit / Length (cm)	31.5
Packed unit / Weight	0.616	Bare Unit / Weight (kg)	0.58
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Colors	Black	Packed unit / Quantity	1

31/05/2026 18:34:29 – For information purposes only, can be subject to change