

# CUFF FOR NEOTENS® ELECTRONIC BLOOD PRESSURE MONITOR - SPENGLER

Ref. 527052



31/05/2026 19:46:06 – For information purposes only, can be subject to change

# CUFF FOR NEOTENS® ELECTRONIC BLOOD PRESSURE MONITOR - SPENGLER

Ref. 527052

Range of cuffs for NEOTENS® electronic blood pressure monitor

## Description

NEOTENS® blood pressure cuffs are very comfortable and come in three sizes:

- Adult cuff M/L (22-42 cm) for NEOTENS® - Reference: 527051
- Child cuff S (16-24 cm) for NEOTENS® - Reference: 527050
- Adult cuff XL (40-55 cm) for NEOTENS® - Reference: 527052

## Features

CE Marking	CE	Medical device	Class I
Bar code	6976054540158	HS Code	90262040
Country	China		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	20	Packed unit / Height (cm)	2
Packed unit / Length (cm)	24	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-25	Max temperature	60
Minimum humidity	NA	Maximum humidity	90

### OTHER FEATURES

31/05/2026 19:46:06 – For information purposes only, can be subject to change