



**Ref. TS53000S**



**CE**  
**0459**



31/05/2026 20:21:48 – For information purposes only, can be subject to change

=====

**Ref. TS53000S**



=====

## Description

=====

## Features

CE Marking	CE0459	Medical device	No
Bar code	3700468215486	HS Code	90189010
Country	France		

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	7	Packed unit / Height (cm)	25
Packed unit / Length (cm)	16	Packed unit / Weight	nc

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	10	Max temperature	70
Minimum humidity	10	Maximum humidity	100

### OTHER FEATURES

 31/05/2026 20:21:48 – For information purposes only, can be subject to change