

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED

®

LUXAMED®

Ref. C1.516.944



+ Graduation

CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



03/06/2026 01:41:25 – For information purposes only, can be subject to change

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED

®

LUXAMED®[Ref. C1.516.944](#)

Innovative pocket-sized dermatoscope with graduated contact plate.

Available in a range of colours.

Description

- Lighting: LED ring with 6 (concentric) LEDs, for a highly detailed image and even lighting.
- Maintenance-free, as no bulb change is required.
- Removable contact plate for easy cleaning
- Continuous zoom with X10 magnification for a distortion-free image
- Optimised head design
- Brightness can be dimmed to avoid glare
- Automatic switch-off after 3 minutes
- LED life: approx. 100,000 hours at constant current
- Illuminance: approx. 10,000 Lux
- Colour temperature: approx. 4,200 K
- Handle made of two materials: aluminium and reinforced fibreglass, compatible with the head of the Luxascope CCT otoscope and ophthalmoscope
- ON/OFF push-button switch
- Power supply: 2x AAA alkaline batteries, included

Available in :

- White: C1.516.914
- Light blue: C1.516.934

 03/06/2026 01:41:25 – For information purposes only, can be subject to change

LUXASCOPE DERMATOSCOPE LED 2.5 V - CONTACT PLATE WITH SCALE - LUXAMED®



Ref. C1.516.944

- Yellow: C1.516.944
- Fuchsia: C1.516.974
- Yellow : C1.516.944
- Kids : C1.916.914
- Black : C1.516.114
- Green : C1.516.954
- Purple : C1.516.964

Features

CE Marking	CE	Medical device	Class I
Bar code	04260381963745	HS Code	90189084
Warranty (months)	24	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Colors	YELLOW		

03/06/2026 01:41:25 – For information purposes only, can be subject to change