

# DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.312

**LUXAMED®**



📅 03/06/2026 02:50:19 – For information purposes only, can be subject to change

## DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.312

**LUXAMED**<sup>®</sup>

Diagnostic Penlight LED OF Luxamed is resistant and lasts over time thanks to LED technology!

This technology provides a shadow-free field of vision and optimum illumination in the area of application. The LEDs are powered by a constant current, guaranteeing considerably longer life.

### Description

Penlight for optimal eye and throat examinations, with LED light approved for eye examinations.

- Top quality ABS plastic housing
- Weight 31.5 g
- Life time LED 100.000 h
- Luminous intensity 18.000 – 22.000 mcd
- Colour temperature: approx. 7000 Kelvin
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

White (Ref : D1.211.912)

Blue (Ref : D1.211.212)

Grey (Ref : D1.211.312)

Yellow (Ref : D1.211.512)

Kids (Ref : D1.311.912)

Black (Ref : D1.211.112)

Red (Ref : D1.211.412)

Green (Ref : D1.211.612)

# DIAGNOSTIC PENLIGHT LED - LUXAMED

Ref. D1.211.312

**LUXAMED®**

## Features

CE Marking	CE	Medical device	Class I
Bar code	04260381961697	HS Code	90189084
Warranty (months)	24	Country	Germany
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.35
Bare Unit / Weight (kg)	0.04		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	10	Maximum humidity	90
<b>OTHER FEATURES</b>			
Colors	Gray		

📅 03/06/2026 02:50:19 – For information purposes only, can be subject to change