

SUNTECH APC BAYONET POCKET CUFFS FOR COMEN NC3/NC5/C30 MONITORS

SUNTECH

Ref. DSU231130U



CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2
2 YEARS
GUARANTEE

03/06/2026 02:49:54 – For information purposes only, can be subject to change

SUNTECH APC BAYONET POCKET CUFFS FOR COMEN NC3/NC5/C30 MONITORS

SUNTECH

Ref. DSU231130U

Suntech APC bayonet armbands are the perfect combination of comfort and reliability!

Description

The special design of the armbands sewn in two pieces offers proven reliability. Designed specifically for all clinical environments, the range of cuffs focuses on ease of application and patient comfort.

- Nylon material with removable pocket
- cuff design based on solid seam extends the service life accompanied by a 2-year warranty.
- Application facilitated by the tapered shape that allows a better wrapping around the arm and thus a better comfort of the patient.

Features

| | | | |
|------------------------------------|------------|---------------------------|---------|
| CE Marking | CE | Medical device | Class I |
| HS Code | 9018905040 | Warranty (months) | 24 |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | nc |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 65 |
| Minimum humidity | 15 | Maximum humidity | 95 |
| OTHER FEATURES | | | |

📅 03/06/2026 02:49:54 – For information purposes only, can be subject to change