

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_B

**ELITE BAGS**



03/06/2026 01:22:06 – For information purposes only, can be subject to change

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_B



**Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner**

### Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening
  
- Material: Polyester gray bitone
- Weight: 0.12 kg.
  
- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03\_N), blue jeans (MEDIAB03\_B) or pink (MEDIAB03\_RC)
  
- Sold empty

# COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_B



## Features

|                                    |               |                         |               |
|------------------------------------|---------------|-------------------------|---------------|
| CE Marking                         | NONE          | Medical device          | No            |
| Bar code                           | 8436564542458 | HS Code                 | 4202.12.99.90 |
| Warranty (months)                  | 24            | Country                 | China         |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |               |
| Packed unit / Width (cm)           | 19            | Bare Unit / Width (cm)  | 6             |
| Packed unit / Height (cm)          | 6             | Bare Unit / Height (cm) | 10            |
| Packed unit / Length (cm)          | 10            | Bare Unit / Length (cm) | 17            |
| Packed unit / Weight               | 0.248         | Bare Unit / Weight (kg) | 0.012         |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |               |
| Minimum temperature                | NA            | Max temperature         | NA            |
| Minimum humidity                   | 30            | Maximum humidity        | 85            |
| <b>OTHER FEATURES</b>              |               |                         |               |
| Colors                             | Bleu jean     | Packed unit / Quantity  | 1             |

03/06/2026 01:22:06 – For information purposes only, can be subject to change