

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_RC

**ELITE BAGS**



03/06/2026 02:14:55 – For information purposes only, can be subject to change

## COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_RC

 **ELITE BAGS**

**Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner**

### Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening
  
- Material: Polyester gray bitone
- Weight: 0.12 kg.
  
- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03\_N), blue jeans (MEDIAB03\_B) or pink (MEDIAB03\_RC)
  
- Sold empty

 03/06/2026 02:14:55 – For information purposes only, can be subject to change

# COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03\_RC



## Features

CE Marking	NONE	Medical device	No
Bar code	8436564542465	HS Code	4202.12.99.90
Warranty (months)	24	Country	China

### DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	19	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	17
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.012

### SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	30	Maximum humidity	85

### OTHER FEATURES

Colors	Pink	Packed unit / Quantity	1
--------	------	------------------------	---

03/06/2026 02:14:55 – For information purposes only, can be subject to change