

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03

ELITE BAGS



03/06/2026 02:15:52 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03



Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner

Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancing device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening

- Material: Polyester gray bitone
- Weight: 0.12 kg.

- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03_N), blue jeans (MEDIAB03_B) or pink (MEDIAB03_RC)

- Sold empty

 03/06/2026 02:15:52 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03



Features

CE Marking	NONE	Medical device	No
Bar code	8436564542199	HS Code	4202129990
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	10
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	6
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	17
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.012
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	30	Maximum humidity	85
OTHER FEATURES			
Colors	Gray	Packed unit / Quantity	1

03/06/2026 02:15:52 – For information purposes only, can be subject to change