

PATIENT CABLE FOR ECG CARDIOMATE * SPENGLER

Ref. 260020



07/06/2026 16:26:20 – For information purposes only, can be subject to change

PATIENT CABLE FOR ECG CARDIOMATE * SPENGLER



Ref. 260020

Patient Cable for ECG CARDIOMATE * SPENGLER

Description

Patient Cable for ECG CARDIOMATE * SPENGLER

Features

CE Marking	NONE	Medical device	No
Bar code	3700446025502	HS Code	8544200090
Warranty (months)	60	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.4

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

📅 07/06/2026 16:26:20 – For information purposes only, can be subject to change