


POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104



📅 07/06/2026 17:22:15 – For information purposes only, can be subject to change

POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104

Those cuffs are autoclavable!

Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

Dimensions of the bladder in cm

Cuff infant (GBS101)	19 x 4 cm
Cuff child (GBS102)	26.5 x 5 cm
Cuff small adult (GBS101)	35.5 x 6 cm
Cuff adult (GBS103)	53 x 6 cm

POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS104

Features

CE Marking	CE	Medical device	Class I
Bar code	3700468203506	HS Code	90189084
Warranty (months)	24	Country	France

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.5

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-10	Max temperature	40
Minimum humidity	30	Maximum humidity	40

OTHER FEATURES

07/06/2026 17:22:15 – For information purposes only, can be subject to change