

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)


Ref. 512200



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



07/06/2026 17:03:26 – For information purposes only, can be subject to change

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Ref. 512200



Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflation bulb. Supplied in a nylon pouch.

Features

| | | | |
|------------------------------------|---------------|---------------------------|-------------------------------|
| CE Marking | CE0459 | Medical device | Class Im (measuring function) |
| Bar code | 3700446007263 | HS Code | 90189010 |
| Warranty (months) | 24 | Country | France |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 14 | Packed unit / Height (cm) | 8 |
| Packed unit / Length (cm) | 20 | Packed unit / Weight | 0.8 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 70 |
| Minimum humidity | 20 | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

📅 07/06/2026 17:03:26 – For information purposes only, can be subject to change