

**VAQUEZ LAUBRY NANO / VELCRO CUFF /
BLACK / CHILD (S)**

Ref. 512210



CE
0459

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
2
ANS
2 YEARS
GUARANTEE

Description

Lightweight and robust Cuff-mounted, it meets the strictest requirements. Large, lightweight, ABS shock-resistant dial, in accordance with standard ISO 81060-1. Non-deformable Oxford nylon cuff with split leather manometer holder. Large inflatio bulb. Supplied in a nylon pouch.

07/06/2026 17:24:34 – For information purposes only, can be subject to change

VAQUEZ LAUBRY NANO / VELCRO CUFF / BLACK / CHILD (S)



Ref. 512210

Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446008628	HS Code	90189010
Warranty (months)	24	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	20	Packed unit / Weight	0.654
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

📅 07/06/2026 17:24:34 – For information purposes only, can be subject to change