

# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE

Ref. 518011



**CE**  
**0459**

DISPOSITIF MEDICAL  
CLASSE  
**Im**  
MEDICAL DEVICE

GARANTIE  
5 ANS  
**5**  
5 YEARS  
GUARANTEE

📅 07/06/2026 17:25:04 – For information purposes only, can be subject to change

## VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

=====

### Description

Modernity and tradition combined. Almost an icon. The most famous of the Vaquez-Laubry®. Available with strap or velcro fastener, in 3 colours and 3 sizes. Poly Cotton cuff with fullgrain leather manometer holder. Large inflation bulb. Supplied in nylon pouch.

Adult (M), cotton, velcro, navy: 518012

Adult (M), cotton, velcro, red: 518013

Adult (M), navy, with straps: 518112

Adult (M), nylon, velcro, grey: 518011

Adult (M), red with straps: 518113

Adult (L), nylon, velcro, grey: 518002

Child/small adult (S), cotton, velcro, motifs : 518022

 07/06/2026 17:25:04 – For information purposes only, can be subject to change

# VAQUEZ-LAUBRY CLASSIC/VELCRO CUFF/LARGE ADULT (L) / NAVY BLUE



Ref. 518011

## Features

CE Marking	CE0459	Medical device	Class Im (measuring function)
Bar code	3700446007287	HS Code	90189010
Warranty (months)	60	Country	France
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	20	Packed unit / Weight	0.75
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	70
Minimum humidity	20	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 07/06/2026 17:25:04 – For information purposes only, can be subject to change